

EEG Toolbox

Developed by Davide Masciola | v1.1-beta

Upload EEG file



Drag and drop file here
Limit 200MB per file • EDF

Browse files

 S001R10.edf 2.5MB ✕

--- DATA EXTRACTER ---

Sampling frequency: 160.0 Hz

Sampling time: 125.0 s

Channels: 64


Channel names: ['Fc5.', 'Fc3.', 'Fc1.', 'Fc2.', 'Fc4.', 'Fc6.', 'C5..', 'C3.

Please select one channel to display (1 to 64):

1−+

Load Channel

Original Signal



Signal Processing

Please select an option

☒ Current Artifact Remover


Select Continent:

- ☒ America (60 Hz)
- ☐ Europe (50 Hz)

Q-factor:

600

601000



☒ Filter

Select Filter:

- ☐ Highpass Filter
- ☐ Lowpass Filter
- ☒ Bandpass Filter
- ☐ Notch Filter

Window Method:

- ☒ Hanning
- ☐ Hamming
- ☐ Blackman
- ☐ Kaiser

Select an EEG band:

- ☐ δ (0.5-3.5 Hz)
- ☒ θ (4-7 Hz)
- ☐ α (8-13 Hz)
- ☐ β (13-30 Hz)
- ☐ γ (>30 Hz)
- ☐ custom

Transition Width Window Method (Hz):




Transition Width Equiripple (Hz):




Transition Width Least Squares (Hz):







Window order:



Equiripple order:



Least Squares order:



Select a FIR Filter:

- ☒ Window Method
- ☐ Equiripple
- ☐ Least Squares

Filtered Signal

